# Pharma Temperature Sensor G3/8"



#### Application/Specified usage

- Aseptic temperature measurement, inline, high precise and fast without product contact
- Temperature measuring in pipes and vessels without opening the process with prefabricated thermowells and build-in systems
- Demounting the sensor without opening the process and without electrical disconnection > avoiding downtime of the equipment at calibration and maintenance!
- Suitable at small pipe diameters with build-in system ESP (available for pipes DN25...DN100)

#### **Application examples**

- Flexible applicable for nearly every temperature measuring task in pipes and vessels
- Safe temperature measuring in hotsteam- and pressure pipes (enclosed process)
- · Monitoring of CIP-/SIP-cleaning

#### Hygienic design/Process connection

- · Hygienic and easy sterilizable installation by using Negele build-in system ESP
- · CIP-/SIP-cleaning up to 140 °C
- · All wetted materials compliant to FDA
- · Sensor completely made of stainless steel
- · 3-A approval for build-in system ESP-G ≥ DN25, ISO20, G1" and ESP-E available

#### Features

- · Short reaction time, very compact measure point
- · Integrated transmitter (optional)
- · Spring mounted gauge slide at TFP-58P
- · Spring mounted sensor tip at TFP-168P and TFP-188P
- · Weight reduced connecting head: non-sensitive to vibrations, hygienic design
- · Electrical connection via M12-plug
- Material (1.4435), inspection certificate 3.1 in scope of delivery (for all product contacting parts)
- · Quick and easy to install with an orbital welding machine
- Temperature sensors and build-in system with predefined and concerted standard lengths reducing product variants and saving storage costs and simplify maintenance
- · Protection class IP 69 K

#### **Options/Accessories**

- · 2 x Pt100 (not retrofittable)
- · 2 x Pt100 with two transmitters (not retrofittable)
- · Programmable transmitter MPU-4 and MPU-M with output 4...20 mA, 2-wire
- · Transmitter Profibus PA and HART protocol
- · Programming adapter MPU-P 9701
- · Transmitter MPU-LCD with integrated display in connecting head
- · Pt100-chip with other classes of accuracy, (1/3 B, 1/10 B)
- · Preassebled cable for M12-plug
- · Fixed cable for TFP-188P in other length and material available

#### **Authorisations**



#### Temperature sensor TFP-58P



Temperature sesor TFP-168P with build-in system ESP-G





Temperature sensor			
Process connection	build-in system ESP	with G3/8" external thread and thermowell	
Insertion length	standard	37 mm, 83 mm, 97 mm, 160 mm	
Materialis connection head		stainless steel 1.4305 (303)	
	protection tube	stainless steel 1.4404	
	cap nut	stainless steel 1.4571	
	spacer	stainless steel 1.4301, Ø 10 mm	
Temperature ranges	ambient	-50+80 °C	
	sensor tip	-50+250 °C	
Operating pressure		50 bar maximum	
Sensing resistor	acc. to DIN EN 60751	1 x Pt100 class A	
Electrical connection	cable gland	M16 x 1.5	
	cable connection	M12-plug 1.4305, 4-pin	
	fixed cable (2.5 m)	PTFE 4 x 0.14 mm2	
Protection type		IP 69 K (with electrical connection M12-plug)	

Transmitter MPU-4, MPU-10, MPU-H, MPU-M					
Temperature ranges	ambient storage	-40+85 °C -55+90 °C			
Measuring ranges	MPU-4, MPU-H, MPU-M	standard: -1040 °C, 050 / 100 / 150 / 200 °C special ranges free programable standard: -200850 °C configuration occurs with Profibus			
Accuracy	input	< ±0.25 °C			
Temperature drift	zero, span	< 0.01 % / K			
Supply	MPU-M, MPU-4 MPU-10 accuracy	835 V DC 932 V DC 0.01 % / V (reference: 12 V DC)			
Output	signal accuracy burden	analog 420 mA (not for MPU-10) < $\pm 0.1$ % of measurement range < $600 \Omega$ (at U <sub>B</sub> = 24 V)			
Humidity	without condensation	098 %			

Accuracy classes of temperature sensors   Tolerances for Pt100 acc. to DIN EN 60751					
Pt100	A	1/3 B	1/10 B		
0°C/100Ω	±0,15 K / ±0,06 Ω	±0,10 K / ±0,04 Ω	±0,03 K / ±0,01 Ω		
100 °C / 138,5 Ω	±0,35 K / ±0,13 Ω	±0,27 K / ±0,10 Ω	±0,08 K / ±0,03 Ω		

Table reaction time	ESF-G-DIN2-10	
Medium temperature 150 °C	t <sub>50</sub>	4,4 s
Medium temperature 150 °C	t <sub>90</sub>	13,1 s

Reaction time

We recommend to use heat-conductive paste. This can reduce the response time up to 50 %.

**Electrical Connection PHARMA** 

#### **Electrical connection without transmitter**

#### With 1 x M12 plug


3

#### Configuration 1st M12 plug



## With M12 plug

### Configuration M12 plug




1: + supply

**Electrical connection with transmitter** 

2: - supply 4...20 mA 3: not connected


#### With 2 x M12 plug



#### Configuration 2nd M12 plug




#### With cable gland



#### With cable gland


#### Configuration strip terminal



#### Electrical connection with two transmitter (TFP-68P)

#### With 1 x M12-plug (sensor 1 + sensor 2)

# Configuration M12-plug



1: + supply (sensor 1)

2: - supply 4...20 mA (sensor 1) 3: - supply 4...20 mA (sensor 2)

4: + supply (sensor 2)

#### With fixed cable




#### Fixed cable connection with 1 x Pt100

wh ye bn gn standard rd rd wh wh PTFE Pt100

#### With 2 x M12-plug (sensor 1)

### Configuration M12-plug



1: + supply (sensor 1)

2: - supply 4...20 mA (sensor 1)


3: not connected 4: not connected

#### Fixed cable connection with 2 x Pt100 (PTFE)

wh 1st Pt100 rd rd ye 2nd Pt100

#### With 2 x M12-plug (sensor 2)






1: + supply (sensor 2)


2: - supply 4...20 mA (sensor 2)


3: not connected


4: not connected

4










TFP-168P / ... / MPU-M



Advices

#### Mechanical connection/Installation



 To guarantee a definite function use the Negele PHARMadapt ESP system.

#### Conventional usage



- · Not suitable for applications in explosive areas.
- Not suitable for applications in security-relevant equipments (SIL).

#### Transport/Storage

5



- · No outdoor storage
- · Dry and dust free
- · Not exposed to corrosive media
- Protected against solar radiation
- · Avoiding mechanical shock and vibration
- Storage temperature -55...+90 °C
- · Relative humidity maximum 98 %

#### Reshipment



- · Sensors shall be clean and must not be contaminated with dangerous media and / or heat-conductive paste.
- Use suitable transport packaging only to avoid damage of the equipment!

#### Cleaning/Maintenance



 In case of using pressure washers, dont't point nozzle directly to electrical connections!

#### Standards and guidelines



 You have to comply with applicable regulations and directives.

#### Disposal



- This instrument is not subject to the WEEE directive 2002/96/EC and the respective national laws.
- Pass the instrument directly on to a specialised recycling company and do not use the municipal collecting points.

#### Advice to EMC



- The device agrees to following standards: EMC directive 2004/108/EG.
- · You have to guarantee the EMC directives for the entire equipement.

#### Conditions for a measuring point according to 3-A-Standard



- The sensors TFP-58P, -68P, -168P, -188P do not require 3-A certification as they do not come into contact with the product.
- · The corresponding PHARMadapt ESP build-in system is 3-A certified.
- Details on the mounting position, self-draining and the position of the leakage hole can be found in the PHARMadapt ESP product information.

## **Temperature Transmitter MPU-LCD with Display**

#### Application/Specified usage

- · 4...20mA transmitter with LCD for Pt100 temperature sensor
- · For installation in temperature sensor
- · Sensor monitoring

#### **Features**

- · 4-digit display with green backlight
- · Temperature measurement in °C and °F
- · Easy range select by one button
- · Lower costs for wiring because of 2-wire technology

#### Note



See product information "MPU-LCD" for details.

#### Accessories

PVC-cable with M12-connection made of 1.4305, IP 69 K, unshielded M12-PVC / 4-5 m PVC-cable 4-pin, length 5 m PVC-cable 4-pin, length 10 m M12-PVC / 4-25 m PVC-cable 4-pin, length 25 m

PVC-cable with M12-connection, brass nickel-plated, IP 67, shielded M12-PVC / 4G-5 m PVC-cable 4-pin, length 5 m PVC-cable 4-pin, length 10 m PVC-cable 4-pin, length 10 m PVC-cable 4-pin, length 25 m

**Programming adapter** 

MPU-P 9701

Programming adapter for MPU-4, MPU-H and MPU-M

### **PVC-cable with M12-connection**

**Option MPU-LCD** 

(display in the connection head)




#### Programming adapter MPU-P 9701



#### **Build-In systems**



Suitable build-in systems for temperature sensors TFP-58p, -68P, 168P, and -188P you will find in product information **Process Connection PHARMadapt ESP**.



#### Temperature sensor version with 2 x Pt100 TFP-58P.2 (connecting head Ø 55 mm, 2 x Pt100, non-sensitive design to vibrations) TFP-68P (like TFP-58P.2, but with higher connecting head and prepared for 2 x transmitter) Sensor length EL in mm 037 059 083 160 **Accuracy class Pt100** Α 1/3B 1/10B Electrical connection only for TFP-58P.2 PG (cable gland M16x1,5) 2 x PG (2 x cable gland M16x1,5) (2 x M12-plug 1.4305) 2 x M12 **Electrical connection only for TFP-68P** M12 (M12-plug 1.4305) 2 x M12 (2 x M12-plug 1.4305) Continue if TFP-68P is selected! No further options for TFP-58P.2! **1st Transmitter** MPU-4 (programmable) Measuring range 1. MPU -10...40 (measuring range -10...40 °C) 0...50 (measuring range 0...+50 °C) 0...100 (measuring range 0...+100 °C) 0...150 (measuring range 0...+150 °C) 0...200 (measuring range 0...+200 °C) хх...уу (special range) **2nd Transmitter** MPU-4 (programmable) Measuring range 2.MPU -10...40 (range -10...40 °C) 0...50 (range 0...+50 °C) 0...100 (range 0...+100 °C) 0...150 (range 0...+150 °C) 0...200 (range 0...+200 °C) хх...уу (special range) TFP-68P/ 083/ M12/ MPU-4/ 0...100/ MPU-4/ 0...100

